

***Sequential same slide multiplex  
immunofluorescence and H&E staining  
for combined phenotypic and  
morphological characterization of  
formalin-fixed paraffin-embedded  
tissue sections***

Kyla Teplitz, Katir K Patel, Amanda Bares, Mael  
Manesse, Mark Burton, Bonnie Phillips, Kate Lillard,  
Anne Hellebust, Melinda Duplessis, Tad George


**indica labs**

**U** Ultivue

**RARECYTE**


# INTRODUCTION

Ultivue's PD-L1 FixVUE™ panel was used for multiplex immunofluorescence staining of CD8, CD68, PD-L1, and pan-Cytokeratin in formalin-fixed, paraffin-embedded (FFPE) samples from 3 serial sections of human tonsil and primary colon and lung tumor biopsies using the Leica® Biosystems BOND RX autostainer. Stained tissues were imaged in five spectrally distinct fluorescence channels (DAPI, FITC, TRITC, Cy5, Cy7) on the RareCyte CyteFinder® II HT Instrument. Slides were de-coverslipped and stained with H&E, then imaged with brightfield using the CyteFinder® II HT Instrument.



# METHODS


Ultivue's PD-L1 FixVUE™ panel was used for multiplex immunofluorescence staining of CD8, CD68, PD-L1, and pan-Cytokeratin in formalin-fixed, paraffin-embedded (FFPE) samples from 3 serial sections of human tonsil and primary colon and lung tumor biopsies using the Leica® Biosystems BOND RX autostainer. Stained tissues were imaged in five spectrally distinct fluorescence channels (DAPI, FITC, TRITC, Cy5, Cy7) on the RareCyte CyteFinder® II HT Instrument. Slides were de-coverslipped and stained with H&E, then imaged with brightfield using the CyteFinder® II HT Instrument.






## ***METHODS***

To segment the tumor and stroma tissue regions, a HALO AI™ classifier was created for the lung and colon H&E images. Fluorescence images were analyzed using the HALO® Highplex FL module to identify CD8+ cytotoxic T-cells, CD68+ macrophages, CD68+/PD-L1+ immuno-suppressive macrophages, pan-CK+ tumor cells, and pan-CK+/PD-L1+ immune-evading tumor cells within the tumor and stromal regions identified by the H&E stain. As a comparison, a classifier was also trained on the fluorescent CK and DAPI signal.






# SEQUENTIAL STAINING AND ANALYSIS WORKFLOW



# ***TUMOR MARKER IS NOT REQUIRED IN PANEL WHEN H&E USED FOR TISSUE ARCHITECTURE CLASSIFICATION***




Cell counts/region comparing a classifier trained on the H&E image vs the fluorescence CK+DAPI signal are comparable. Both classifiers provide equivalent results, indicating that the tumor marker could be dropped from the fluorescence panel if using the H&E for tissue architecture classification.

\*Lung 6 tumor/stroma segmentation was applied for Lung 5, since the tissue was damaged during deconvolving.




## COLON TISSUE CLASSIFIED FROM H&E IMAGE



(A) H&E image of colon tissue (B) classified into tumor epithelium (red), stroma (green), and normal epithelium (blue) regions. (C) Tumor region and (D) normal region of the fluorescence image shown in closeup and segmented with color-coded phenotypes outlines: CK-positive (cyan), CD8-positive (green), CD68- positive (red), PD-L1-positive (pink), negative cytoplasm (grey).



## ***LUNG TISSUE CLASSIFIED FROM H&E IMAGE***



(A) H&E image of lung tissue (B) classified into tumor epithelium (red), and stroma (green). (C) Tumor region and (D) stromal region of the fluorescence image shown in closeup and segmented with color-coded phenotypes outlines: CK-positive (cyan), CD8-positive (green), CD68-positive (red), PD-L1-positive (pink), negative cytoplasm (grey).



# CONCLUSIONS


Here we demonstrate a tissue-preserving workflow to generate H&E images from a slide that is previous stained and imaged in fluorescence on the CyteFinder® II HT Instrument.

These H&E images can be used to delineate tissue architectural regions such as tumor and stroma, eliminating the need for a tumor-specific biomarker in the fluorescence panel.

H&E trained classifiers perform equivalently to classifiers trained on fluorescent signal from DAPI and pan-CK.

Ultivue's PD-L1 FixVUE™ panel produces reproducible staining results across multiple serial sections.

HALO® Highplex FL can be used to identify phenotypic cell populations and expression levels within tissue architectural regions.



# ***PRODUCT INFORMATION***


**CyteFinder® II HT Instrument**

**PD-L1 FixVUE™ panel**

**HALO® image analysis platform**

**HALO® Highplex FL**

**RARECYTE®** **indica labs** **U** **Ultivue**





## CONTACT INFORMATION

- To learn more about PD-L1 FixVUE™ panel, email [contact@ultivue.com](mailto:contact@ultivue.com)
- To learn more about the CyteFinder® II HT Instrument, email [info@rarecyte.com](mailto:info@rarecyte.com)
- To learn more about HALO® image analysis workflows, email [info@indicalab.com](mailto:info@indicalab.com)

