Abstract #2629







¹Clinical Pharmacodynamics Biomarker Program, Applied/Developmental Research Directorate, Frederick National Laboratory for Cancer Research, Frederick, MD 21702; ²Ultivue, Inc., Cambridge, MA 02138; ³Developmental Therapeutics Clinic/Early Clinical Trials Development Program, Division of Cancer Treatment and Diagnosis, National Cancer Institute, Bethesda, MD 20892; ⁵Center for Cancer Research, National Cancer Institute, Bethesda, MD 20892; ⁴Division of Cancer Institute, Bethesda, MD 20892; ⁵Center for Cancer Research, National Cancer Institute, Bethesda, MD 20892; ⁴Division of Cancer Institute, Bethesda, MD 20892; ⁴Division of Cancer Treatment and Diagnosis, National Cancer Institute, Bethesda, MD 20892; ⁴Division of Cancer Institute, Be

Abstract

<u>Background</u>: Immune checkpoint inhibitors promote antitumor immune responses by enhancing T-cell activity. Measuring the pharmacodynamic effects of these drugs is challenging, as it requires assessing both immune cell and cancer cell populations. To evaluate T-cell activation in tumor tissue from patient biopsies, we developed a robust multiplexed immunofluorescence assay.

Methods: Our assay uses novel oligo-conjugated antibodies (Ultivue) for simultaneous quantitation of TCR activation (phospho-CD3zeta), immune checkpoint signaling via PD-1 (p-SHP1/p-SHP2), and the net stimulation/inhibition resulting from the integration of these two pathways in CD8 cells (p-ZAP70), while also providing the proximity of CD8 cells to tumor tissues, identified by β -catenin. The method was clinically validated using custom tissue microarrays (TMA) containing tumor biopsies of 3 different histologies (CRC, NSCLC, and breast).

<u>Results</u>: From a total of 192 tumor core biopsies, 20/64 NSCLC, 9/64 CRC, and 3/65 breast TMA cores were found to have a significant number of CD8+ tumor infiltrating lymphocytes (TILs) at baseline (>50 cells in the examined section) In 18 of the 20 NSCLC cores, ≥50% of CD8 cells both inside and outside of the tumor were activated (CD3z-pY142+). In 6/9 CRC cores, ≥50% of CD8+ cells inside tumor tissues were activated, and in 4/9 CRC cores, ≥50% of CD8+ cells in stroma were activated. In 2/3 breast tumor cores, 90% of CD8+ cells inside tumor tissues were activated; in the remaining core, 90% of CD8+ cells in stroma were activated. Interestingly, all 192 cores had minimal to no expression of activated Zap70 (pY493) in CD8+ cells

<u>Conclusions</u>: Depending on tumor histology, baseline biopsy samples may contain variable numbers of activated CD8+ TILs (CD3z-pY142+), which may reside inside or outside of tumor regions and express very low levels of Zap70-pY493. Anti-PD-1 therapy is predicted to enhance T-cell cytotoxic activity, as demonstrated by an increased number of TILs and elevated Zap70-pY493 expression. This assay is being used for pharmacodynamic evaluations in ongoing immunotherapy clinical trials. Funded by NCI Contract No HHSN261200800001E.

Background



PD-1 modulation of TCR signaling. Binding of PD-L1 ligands to PD-1 leads to the binding of SHP-2 to phosphorylated ITSM and overall inhibition of T cell receptor (TCR) signaling through blockade of CD3z chain phosphorylation and Zap-70 association.

Materials and Methods

Development of oligo-conjugated 5-plex quantitative immunofluorescent assay (IFA): Validated antibodies to CD8, CD3z-pY142, Zap70-pY493, and β -catenin were conjugated to specific oligonucleotides and detected by complementary fluor-conjugated oligos (FITC, TRITC, Cy5, and Cy7, respectively; Ultivue). DAPI was included in the panel to assess cellularity by nuclear staining. Control tissue and cell pellet slides, and human tumor TMAs (Indivumed), were stained on Leica Bond; images were acquired on a Zeiss Axioscanner and analyzed by Definiens Software. **β**-catenin



Tissue and cell pellet controls for quantitative multiplex IFA (qmIFA) development: β -catenin staining of tumor tissue for tumor segmentation was validated using an MTU951 Multi-tumor Tissue Microarray (US Biomax). CD3z-pY142 and CD8 tissue staining were validated on human tonsil. Zap70-pY493 staining was validated on anti-CD3/CD28 beadactivated T-cell pellets. The 5-plex IFA staining (including DAPI) was clinically validated on human tumor tissue microarrays from 3 different histologies: CRC, NSCLC, and breast (Indivumed).

Quantitative analysis of 5-plex IFA: Image analysis was performed in Definiens Architect 2.4.2 Tissue Studio IF. Tumor tissue was identified by β-catenin staining, and stroma was identified by absence of β-catenin staining. The tissue of each core was classified as tumor or stroma using the Composer Training in ROI detection. In Cellular Analysis, nuclei were detected and assigned to either tumor or stroma. Cell Simulation was based on cytoplasmic staining. In addition, TILs were identified by CD8 or CD3z-pY142 staining; cells that were CD8+ CD3z-pY142+, CD8+ CD3z-pY142-, or CD8-CD3z-pY142+ were identified in Cell Classification using the coexpression feature. These cells were assigned to either the tumor or stroma based on B-catenin staining

Definiens IO-PD Multiplex IFA Solution Overview



9 Prostate





















A multiplex immunofluorescence assay to assess immune checkpoint inhibitor-targeted **CD8** activation and tumor co-localization in FFPE tissues

Tony Navas¹, Kristin Fino¹, King Leung Fung¹, Facundo Cutuli¹, Robert J. Kinders¹, Aditi Sharma², Geraldine O'Sullivan Coyne³, Alice P. Chen³, Toby Hecht⁴, James H. Doroshow^{3,5}, Ralph E. Parchment¹

Results



Figure 5. Definiens Image Analysis and Quantitation of CRC, NSCLC, and **Breast TMA Cores**







Frederick National Laboratory

for Cancer Research

CD8+ CD3z-pY142-CD8+ CD3z-pY142+ CD8- CD3z-pY142+

Table 1. Quantitation of total number of activated CD8 T cells in relation to tumor tissue or surrounding stroma in select CRC and breast TMA cores

CRC C2	Biopsy	Stroma	% Stromal Cells	Tumor	% Tumor Cells
Total Cells	6,783	2,841	42.0%	3,942	58.0%
CD8+ CD3z-	61	31	1.10%	30	0.80%
CD8+ CD3z+	592	43	1.50%	549	13.90%
CD8- CD3z-	5575	2,638	34.9%	2,937	32.5%
CD8- CD3z+	555	129	4.50%	426	10.80%
CRC C7	Biopsy	Stroma	% Stromal Cells	Tumor	% Tumor Cells
Total Cells	5,132	2,762	53.8%	2,370	46.2%
CD8+ CD3z-	2	2	0.07%	0	0.00%
CD8+ CD3z+	183	145	5.25%	38	1.60%
CD8- CD3z-	4587	2,340	38.5%	2,247	41.0%
CD8- CD3z+	360	275	9.96%	85	3.59%
NSCLC D6	Biopsy	Stroma	% Stromal Cells	Tumor	% Tumor Cells
Total Cells	5,330	1,834	34.4%	3,496	65.6%
CD8+ CD3z-	236	61	3.30%	175	5.00%
CD8+ CD3z+	582	103	5.60%	479	13.70%
CD8- CD3z-	4312	1,572	20.2%	2,740	44.0%
CD8- CD3z+	200	98	5.30%	102	2.90%
NSCLC E8	Biopsy	Stroma	% Stromal Cells	Tumor	% Tumor Cells
NSCLC E8 Total Cells	Biopsy 5,868	Stroma 3,432	% Stromal Cells 58.5%	Tumor 2,436	% Tumor Cells 41.5%
NSCLC E8 Total Cells CD8+ CD3z-	Biopsy 5,868 70	Stroma 3,432 40	% Stromal Cells 58.5% 1.20%	Tumor 2,436 30	% Tumor Cells 41.5% 1.20%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+	Biopsy 5,868 70 636	Stroma 3,432 40 472	% Stromal Cells 58.5% 1.20% 13.80%	Tumor 2,436 30 164	% Tumor Cells 41.5% 1.20% 6.70%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z-	Biopsy 5,868 70 636 4278	Stroma 3,432 40 472 2,229	% Stromal Cells 58.5% 1.20% 13.80% 23.4%	Tumor 2,436 30 164 2,049	% Tumor Cells 41.5% 1.20% 6.70% 25.7%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+	Biopsy 5,868 70 636 4278 884	Stroma 3,432 40 472 2,229 691	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10%	Tumor 2,436 30 164 2,049 193	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+ Breast C6	Biopsy 5,868 70 636 4278 884 Biopsy	Stroma 3,432 40 472 2,229 691 Stroma	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells	Tumor 2,436 30 164 2,049 193 Tumor	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+ Breast C6 Total Cells	Biopsy 5,868 70 636 4278 884 Biopsy 5,966	Stroma 3,432 40 472 2,229 691 Stroma 3,162	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53%	Tumor 2,436 30 164 2,049 193 Tumor 2,804	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z-	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53% 1.08%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z+	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z-	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z- CD8+ CD3z+ CD8- CD3z- CD8- CD3z+	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast E1	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311 Biopsy	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261 Stroma	% Stromal Cells 58.5% 1.20% 13.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3% % Stromal Cells	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50 Tumor	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8% % Tumor Cells
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ CD8- CD3z+ Breast E1 Total Cells	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311 Biopsy 8,058	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261 Stroma 2,237	% Stromal Cells 58.5% 1.20% 1.3.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3% % Stromal Cells	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50 Tumor 5,821	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8% % Tumor Cells 72.2%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast E1 Total Cells CD8+ CD3z-	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311 Biopsy 8,058 66	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261 Stroma 2,237 29	% Stromal Cells 58.5% 1.20% 1.3.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3% % Stromal Cells 28% 1.30%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50 Tumor 5,821 37	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8% % Tumor Cells 72.2% 0.64%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast E1 Total Cells CD8+ CD3z- CD8+ CD3z-	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311 Biopsy 8,058 66 333	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261 Stroma 2,237 29 128	% Stromal Cells 58.5% 1.20% 1.3.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3% % Stromal Cells 1.08% 1.08% 3.64% 40.0% 8.3% % Stromal Cells 28% 1.30% 5.72%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50 Tumor 5,821 37 205	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8% % Tumor Cells 72.2% 0.64% 3.52%
NSCLC E8 Total Cells CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast C6 Total Cells CD8+ CD3z- CD8+ CD3z- CD8+ CD3z+ CD8- CD3z+ Breast E1 Total Cells CD8+ CD3z- CD8+ CD3z- CD8+ CD3z- CD8+ CD3z-	Biopsy 5,868 70 636 4278 884 Biopsy 5,966 46 150 5459 311 Biopsy 8,058 66 333 7140	Stroma 3,432 40 472 2,229 691 Stroma 3,162 34 115 2,752 261 Stroma 2,237 29 128 1,802	% Stromal Cells 58.5% 1.20% 1.3.80% 23.4% 20.10% % Stromal Cells 53% 1.08% 3.64% 40.0% 8.3% % Stromal Cells 28% 1.30% 5.72% 8.4%	Tumor 2,436 30 164 2,049 193 Tumor 2,804 12 35 2,707 50 Tumor 5,821 37 205 5,338	% Tumor Cells 41.5% 1.20% 6.70% 25.7% 7.90% % Tumor Cells 47.0% 0.43% 1.25% 43.5% 1.8% % Tumor Cells 3.5% 1.8% % Tumor Cells 63.9%

Table 2. Summary of T cells found in CRC, NSCLC, and breast TMAs

	Number of cores	TMA cores with ≥50 total CD8+ Tcells	Stage	TMA cores with ≥50% of the total activated T cells in location	
Histology				Tumor	Stroma
Breast	65	3	-	2	1
Colorectal	64	4	II	2	2
		5		4	1
NSCLC	64	20	II-IV	17	3

Summary and Clinical Implications

- We have developed a robust quantitative multiplex IO-PD immunofluorescence assay that quantitatively detects CD8+ cells and their activation status in relation to tumor tissues, as delineated by β -catenin. To analyze stained tissue, we have also developed algorithms by Definiens Architect to quantify activated CD8+ T cells (CD3z-pY142+) both inside tumor tissues as well as in surrounding stroma using CRC, NSCLC, and breast TMA samples from Indivumed.
- Based on the TMA tumor staining, 3/65 breast TMA cores had >50 total infiltrating lymphocytes in the biopsy, of which 2 cores had 90% of activated CD8+ T cells inside the tumor and 1 core had 90% of activated CD8+ T cells outside the tumor. There was minimal Zap70-pY493 expression in T cells.
- In the CRC TMA, 8/64 cores had >50 total infiltrating lymphocytes in the biopsy, of which 6/8 cores had >50% of activated CD8+ T cells inside the tumor and 2/8 cores had >80% of activated CD8+ T cells outside the tumor. Again, no Zap70-pY493 expression in T cells was found at baseline.
- In the NSCLC TMA, 20/64 cores had >50 total infiltrating lymphocytes in the biopsy, of which 17/20 cores had >50% of activated CD8+ T cells inside the tumor and 18/20 cores had >50% of activated CD8+ T cells outside the tumor. As in the other histologies, there was no Zap70-pY493 expression detected in T cells.
- Depending on tumor histology, baseline biopsy samples may contain variable numbers of activated CD8+ TILs (CD3zpY142+), which may reside inside or outside of tumor regions and express very low levels of Zap70-pY493. This assay is being used for pharmacodynamic evaluations in ongoing immunotherapy clinical trials. The assay will be made available to the public via https://dctd.cancer.gov/ResearchResources-biomarkers.htm.

Acknowledgements

Acknowledgments: We thank Jamie Buell, Sean Downing, and Louis Levy for their assistance with oligo-conjugated multiplex IF assay panel development (Ultivue). We thank Brad Gouker for histotechnological services. We thank Rachel Andrews, Manisha Mohandoss and Gabe Benton, Leidos Biomed QC Group, for providing qualified materials and reagents during assay development. We thank Dr. Laura K. Fogli, Kelly Services, for editorial assistance in the preparation of this poster.



Funding: This project has been funded in whole or in part with federal funds from the National Cancer Institute, National Institutes of Health, under Contracts No. HHSN261200800001E. The content of this publication does not necessarily reflect the views or policies of the Department of Health and Human Services, nor does mention of trade names, commercial products, or organizations imply endorsement by the U.S. Government.